

PIP4K2A Polyclonal Antibody

Catalog number: 12469-1-AP

Size: 20 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PIP4K2A; FLJ13267; PI5P4KA;

PIP5K2A; PIP5KII-alpha;

PIP5KIIA; PIPK



mouse brain tissue were subjected to SDS PAGE followed by western blot with 12469-1-AP(PIP4K2A Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse brain tissue
Calculated PIP4K2A MW:	406aa,46kd
Observed PIP4K2A MW:	46kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3152
GenBank accession number:	BC018034
Gene ID (NCBI):	5305
Full name:	Phosphatidylinositol-5-phosphate 4-kinase, type II, alpha

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.