

BL21(plysS) -lysate Polyclonal Antibody

Catalog number: 10562-1-Ig

Size: 105 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

plyss lysate

Applications

Tested applications:	ELISA, WB
Species specificity:	Bacterial; other species not tested.
Calculated BL21(plysS) - lysate MW:	1-200 kDa
Observed BL21(plysS) -lysate MW:	
Positive WB detected in	Total bacteria lysates of BL21(plysS)
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag0045
GenBank accession number:	
Gene ID (NCBI):	
Full name:	Undefined

Product information

Purification method:	Protein A purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.