

Save over \$1,000 on the next generation of Wiley Mini-Mills!



Single Speed



Variable Speed

- Power: 115V 60Hz or 230V 50Hz
- Speed: 1200 rpm free spindle speed
- Chamber Size: 1.575" Diameter
- Stationary Blades: 2
- Rotating Blades: 4
- Max Sample Input Size: .5" Diameter
- Continuous Torque: 1.4 Nm
- Exterior Material: Powder Coated Steel

- Power: 115V 60Hz or 230V 50Hz
- Speed: 1-2400 rpm free spindle speed
- Chamber Size: 1.575" Diameter
- Stationary Blades: 2
- Rotating Blades: 4
- Max Sample Input Size: .5" Diameter
- Continuous Torque: 1.4 Nm
- Exterior Material: Powder Coated Steel

The Mini Mills are all steel Benchtop Cutting Mills featuring a 1.575-inch diameter sealed cutting chamber with two adjustable stationary blades and a rotating cutting head with four cutting edges that revolve at high speed to produce a shearing action that minimizes heat build-up.

The hard, steel cutting edges on the knives allow milling of a wider range of materials, including plastics, animal feed, bones, cables, cardboard, electronic scrap, feed pellets, foils, food, leather, material mixtures, non-ferrous metals, paper, pharmaceutical products, refuse-derived fuels, rubber, waste, and wood.

Thomas No.	Description	List Price	Sale Price
23A00D287	Wiley® Mini Mill, Single-Speed, 115V 60HZ	\$6,998.96	\$5,988.73
23A00D288	Wiley® Mini Mill, Single-Speed, 230V 50HZ	\$8,142.61	\$7,086.67
23A00D289	Wiley® Mini Mill, Variable-Speed, 115V 60HZ	\$10,883.64	\$9,557.75
23A00D290	Wiley® Mini Mill, Variable-Speed, 230V 50HZ	\$11,516.91	\$10,114.86

Use
PROMO Code:
0324
at checkout

Special Pricing Valid Until June 30, 2024!