

Dishes for general-purpose weighing, dispensing, storage, food processing, evaporating or drying

## Aluminum Weighing Dish

- Designed for moisture determination balances with slightly tapered walls for stacking



## Oil Free Aluminum Weighing Boats


**NEW**

Passed through a rigorous, three-step cleansing process to remove forming oil or residue for more reliable weighing and moisture sample analysis.

- Processed and packaged in clean room environment
- Lot traceable



Aluminum Weighing Dish



120934  
120934FR  
15 mL



120938  
120938FR  
45 mL



HS14522  
HEA14522FR  
40 mL



120933  
120933FR  
1 mL



120940  
120940FR  
0.25 mL

## Aluminum Weighing Dish

Item No.	Volume mL	Dimensions in	Dimensions mm	Packaging	Tab	Oil Free
120938	45 mL	3.5 x 0.3	90 x 7	50		
120938FR	45 mL	3.5 x 0.3	90 x 7	50		x
HS14522	40 mL	4 x 0.3	101 x 8	50	x	
HEA14522FR	40 mL	4 x 0.3	101 x 8	50	x	x
120934	15 mL	2.7 x 0.2	70 x 6	50		
120934FR	15 mL	2.7 x 0.2	70 x 6	50		x
120933	1 mL	0.98 x 0.3	25 x 7	50		
120933FR	1 mL	0.98 x 0.3	25 x 7	50		x
120940	0.25 mL	0.5 x 0.1	13 x 3.5	100	x	
120940FR	0.25 mL	0.5 x 0.1	13 x 3.5	100	x	x