

## PLS3 Polyclonal Antibody

Catalog number: 12917-1-AP


Size: 51 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PLS3; Plastin 3, plastin 3 (T isoform), PLS3, T plastin



A431 cells were subjected to SDS PAGE followed by western blot with 12917-1-AP ( PLS3 Antibody) at dilution of 1:600

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PLS3 MW:	630aa,71 kDa
Observed PLS3 MW:	68 kDa
Positive WB detected in	A431 cells, Transfected HEK-293 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3574
GenBank accession number:	BC039049
Gene ID (NCBI):	5358
Full name:	Plastin 3 (T isoform)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.