

PPIL4 Polyclonal Antibody

Catalog number: 12538-1-AP


Size: 39 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

PPIL4; HDCME13P



HeLa cells were subjected to SDS PAGE followed by western blot with 12538-1-AP (PPIL4 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	HeLa cells, HEK-293 cells, mouse kidney tissue
Calculated PPIL4 MW:	492aa, 57kd
Observed PPIL4 MW:	58-65kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3249
GenBank accession number:	BC016984
Gene ID (NCBI):	85313
Full name:	Peptidylprolyl isomerase (cyclophilin)-like 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.