

Specification Sheet

PRODUCT: 10X Tris Buffered Saline

CATALOG NO.: 472845 and
492845

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 1 Year from the date
of manufacture

STORAGE CONDITIONS: 4 - 40°C

Test Results

| Test | Specification | Test Method | Result |
|-------------------------------|------------------------|----------------------------|--------|
| Appearance | Clear/Colorless Liquid | Visual | TBD |
| pH | 7.6 ± 0.1 | Potentiometric | TBD |
| Osmolarity (1:10 Dilution) | 280 ± 20 mOsm | Freeze Point Depression | TBD |
| Protease | None Detected | Enzymatic | TBD |
| DNase | None Detected | Fluorometric | TBD |
| RNase | None Detected | Fluorometric | TBD |
| Peroxide | ≤ 5ppb | Enzymatic | TBD |
| Sterility | No Growth | USP <71> | TBD |
| Endotoxins | ≤ 1 EU/mL | USP <85> | TBD |

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538

Simple Lifecycle