

## DISP1 Polyclonal Antibody

Catalog number: 12041-1-AP

Size: 54 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

DISP1; DISP1, DISPA,

DKFZp434I0428, FLJ43740,

Protein dispatched homolog 1



HeLa cells were subjected to SDS PAGE followed by western blot with 12041-1-AP (DISP1 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human; other species not tested.
Calculated DISP1 MW:	1524aa, 171 kDa
Observed DISP1 MW:	150-171 kDa
Positive WB detected in:	HeLa cells, HEK-293 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag2670
GenBank accession number:	BC011542
Gene ID (NCBI):	84976
Full name:	Dispatched homolog 1 (Drosophila)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.