



Medical Procedure Mask with Earloops

safeForce® procedure masks combine high breathability with exceptional filtration efficiency, while offering fluid resistance against bodily fluids and others. This 3-layer pleated mask is hypoallergenic and features an adjustable nose piece and non-latex earloop for easy donning and added comfort.

Meets ASTM Level 2 Requirements:

- BFE (Bacterial Filtration Efficiency): **>99.9%**
- PFE (Particle Filtration Efficiency): **>98%**
- Fluid Resistance: **120 mm Hg**
- Differential Pressure: **4.5 mm H₂O/cm²**
- Flammability: **Class 1**
- FDA 510(K) Registered

Packaging:

50 masks/Box; 40 Boxes/Case; 2000 Ea/Case



Thomas No.	Description	Quantity
20A00K947	Medical Procedure Mask with Earloops, 50 masks/Box; 40 Boxes/Case; 2000 Ea/Case	50/Pk 2000/Cs

Disclaimer: Pricing and offers subject to change without notice.

Connect With Us:



ThomasSci.com

General Lab 833.544.7447 ■ Production 833.750.1234

Thomas
Scientific