

PLEKHA5 Polyclonal Antibody

Catalog number: 13045-1-AP

Size: 71 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

PLEKHA5; FLJ10667; FLJ31492;

KIAA1686; PEPP2



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13045-1-AP (PLEKHA5 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse testis tissue, human heart tissue, mouse kidney tissue
Calculated PLEKHA5 MW:	1116aa, 127kd
Observed PLEKHA5 MW:	127kd, 74kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3575
GenBank accession number:	BC013133
Gene ID (NCBI):	54477
Full name:	Pleckstrin homology domain containing, family A member 5

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.