

ZNF443 Polyclonal Antibody

Catalog number: 13667-1-AP

Size: 20 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

ZNF443; zinc finger protein 443,

ZK1, ZNF443



mouse ovary tissue were subjected to SDS PAGE followed by western blot with 13667-1-AP(ZNF443 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse; other species not tested.
Calculated ZNF443 MW:	671aa,78 kDa
Observed ZNF443 MW:	72kd
Positive WB detected in	Mouse ovary tissue, mouse lung tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4585
GenBank accession number:	BC032753
Gene ID (NCBI):	10224
Full name:	Zinc finger protein 443

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.