

TTC4 Polyclonal Antibody

Catalog number: 11878-1-AP

Size: 49 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:


TTC4; DKFZp781B0622,

FLJ41930, TPR repeat protein 4,

TTC4



Immunohistochemical of paraffin-embedded human lymphoma using 11878-1-AP(TTC4 antibody) at dilution of 1:50 (under 10x lens)



HeLa cells were subjected to SDS PAGE followed by western blot with 11878-1-AP(TTC4 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated TTC4 MW:	387aa, 45 kDa
Observed TTC4 MW:	45 kDa
Positive WB detected in:	HeLa cells
Positive IHC detected in:	Human lymphoma tissue
Recommended dilution:	WB: 1:500-1:5000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2474
GenBank accession number:	BC001276
Gene ID (NCBI):	7268
Full name:	Tetratricopeptide repeat domain 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.