

FBXO18 Polyclonal Antibody

Catalog number: 14275-1-AP

Size: 21 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

FBXO18; FBH1; FLJ14590;

Fbx18; MGC131916;

MGC141935; MGC141937



PC-3 cells were subjected to SDS PAGE followed by western blot with 14275-1-AP (FBXO18 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	PC-3 cells, HeLa cells
Calculated FBXO18 MW:	118kd
Observed FBXO18 MW:	118-123kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag5622
GenBank accession number:	BC032674
Gene ID (NCBI):	84893
Full name:	F-box protein, helicase, 18

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.