



## Rogosa SL Broth

M407

Rogosa SL Broth is used as a selective medium for cultivation of oral, vaginal and faecal Lactobacilli.

### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	10.000
Yeast extract	5.000
Dextrose	10.000
Arabinose	5.000
Saccharose	5.000
Sodium acetate	15.000
Ammonium citrate	2.000
Monopotassium phosphate	6.000
Magnesium sulphate	0.570
Manganese sulphate	0.120
Ferrous sulphate	0.030
Polysorbate 80	1.000
Final pH ( at 25°C)	5.4±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 59.72 grams in 1000 ml distilled water. Heat, if necessary, to dissolve the medium completely. Add 1.32 ml glacial acetic acid. Mix thoroughly, distribute into culture tubes or flasks. Heat to 90 - 100°C for 2-3 minutes. Cool to 45°C for direct inoculation. DO NOT AUTOCLAVE.

### Principle And Interpretation

Rogosa SL Broth, is known as RMW Broth, it is a modification of media formulated by Rogosa, Mitchell and Wiseman (3, 4). This media is used for isolation, enumeration and identification of Lactobacilli from foodstuffs and clinical specimens (1, 2). Accompanying bacterial flora is suppressed due to the low pH of the medium and also because of the high sodium acetate concentration.

Casein enzymic hydrolysate, yeast extract provide nitrogenous compounds, sulphur, trace elements and vitamin B complex, essential for growth of Lactobacilli. Dextrose, Arabinose, Saccharose are the fermentable carbohydrates. Polysorbate 80 is the source of fatty acids. Ammonium citrate and sodium acetate inhibit moulds, Streptococci and many other organisms. Monopotassium phosphate provides buffering capability. Magnesium sulphate, manganese sulphate and ferrous sulphate are sources of inorganic ions. Low pH of the medium and addition of acetic acid makes the medium selective for Lactobacilli inhibiting other bacterial flora (2).

It is recommended that the plates should be incubated at 30°C for 5 days or at 37°C for 3 days in an atmosphere of 95% hydrogen and 5% carbon dioxide (5). High acetate concentration and acidic pH suppress many strains of other lactic acid bacteria. The salt in the formulation makes the medium unsuitable for isolation of dairy lactobacilli e.g. *L. lactis*, *L. bulgaricus* and *L. helveticus* (2, 4).

### Quality Control

#### Appearance

Cream to yellow homogeneous soft lumps which can be easily broken down to powder form.

#### Colour and Clarity of prepared medium

Light yellow coloured clear to slightly opalescent solution in tubes

#### Reaction

Reaction of 6.0% w/v aqueous solution with 0.132% acetic acid at 25°C pH : 5.4±0.2

**pH**

5.20-5.60

**Cultural Response**

M407: Cultural characteristics observed after an incubation at 35-37°C for 40-48 hours, in presence of 5% Carbon dioxide(CO<sub>2</sub>) and 95% H<sub>2</sub>.

<b>Organism</b>	<b>Inoculum (CFU)</b>	<b>Growth</b>
<i>Lactobacillus casei</i> ATCC 9595	50-100	good - luxuriant
<i>Lactobacillus fermentum</i> ATCC 9338	50-100	good - luxuriant
<i>Lactobacillus leichmanni</i> ATCC 4797	50-100	good - luxuriant
<i>Lactobacillus plantarum</i> ATCC 8014	50-100	good - luxuriant
<i>Staphylococcus aureus</i> ATCC 25923	≥10 <sup>3</sup>	inhibited

**Storage and Shelf Life**

Store dehydrated powder at 2-8°C in tightly closed container and use freshly prepared medium. Use before expiry date on the label.

**Reference**

- Downes F. P. and Ito K., (Eds.), Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
- MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification- Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore. Md.
- Rogosa M., Mitchell J. A. and Wiseman R. F., 1951, J. Bacteriol., 62, 132-133.
- Rogosa M., Mitchell J. A. and Wiseman R. F., 1951, J. Dental Res. 30:682.
- Sharpe M. L. (Ed.), 1960, Lab-Practice, 9(4): 223.

Revision : 1 / 2011

**Disclaimer :**

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.