

## DBT Polyclonal Antibody

Catalog number: 12451-1-AP

Size: 37 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

DBT; BCATE2; E2; E2B;

MGC9061



human brain tissue were subjected to SDS PAGE followed by western blot with 12451-1-AP(DBT Antibody) at dilution of 1:3000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human brain tissue, A375 cells, K-562 cells
Calculated DBT MW:	482aa, 53kd
Observed DBT MW:	53kd
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP =

Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3101
GenBank accession number:	BC016675
Gene ID (NCBI):	1629
Full name:	Dihydrolipoamide branched chain transacylase E2

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.