

SELK Polyclonal Antibody

Catalog number: 12561-1-AP

Size: 23 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

SELK; HSPC030; HSPC297;

MGC17057



K-562 cells were subjected to SDS PAGE followed by western blot with 12561-1-AP (SELK Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	K-562 cells
Calculated SELK MW:	94aa, 10kd
Observed SELK MW:	10kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3143
GenBank accession number:	BC013162
Gene ID (NCBI):	58515
Full name:	Selenoprotein K

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.