

## MAFF Polyclonal Antibody

Catalog number: 12771-1-AP

Size: 63 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:


MAFF; hMafF, MAFF,

Transcription factor MafF, U

MAF



Immunohistochemical of paraffin-embedded human ovary tumor using 12771-1-AP(MAFF antibody) at dilution of 1:100 (under 10x lens)



human lung tissue were subjected to SDS PAGE followed by western blot with 12771-1-AP(MAFF antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated MAFF MW:	17.8 kDa
Observed MAFF MW:	18-23 kDa
Positive WB detected in:	Human lung tissue
Positive IHC detected in:	Human ovary tumor tissue
Recommended dilution:	WB: 1:500-1:5000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3460
GenBank accession number:	BC015037
Gene ID (NCBI):	23764
Full name:	V-maf musculoaponeurotic fibrosarcoma oncogene homolog F (avian)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.