

PCTAIRE3 Polyclonal Antibody

Catalog number: 12207-1-AP

Size: 46 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PCTK3; CDK18, PCTAIRE,
PCTAIRE motif protein kinase 3,

PCTAIRE protein kinase 3,

PCTAIRE3, PCTK3



human brain tissue were subjected to SDS PAGE followed by western blot with 12207-1-AP(PCTK3 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PCTAIRE3 MW:	472aa,54 kDa
Observed PCTAIRE3 MW:	54 kDa
Positive WB detected in	Human brain tissue, HEK-293 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag2851
GenBank accession number:	BC011526
Gene ID (NCBI):	5129
Full name:	PCTAIRE protein kinase 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.