



### Specifications

Electrical requirements	110VAC 60Hz single phase, 900 watts
Usable chamber dimensions	16.5 x 13.25 x 13.75" (WxDxH)
Unit dimensions	28 x 24 x 27" (WxDxH)
Shipping dimensions	32 x 28 x 27" (WxDxH)
Temperature range	Ambient temperature to 400°F
Temperature uniformity	+/- 2°F
Temperature controller precision	0.1°F
Timer range	1 to 9999 minutes
Racks included	2 (additional racks available at extra cost)
Maximum rack capacity	5 (height between each rack is 2 inches)
Exhaust outlet diameter	1 1/4"
Unit weight	95 Lb
Shipping weight	150 Lb
Access port with cover	Optional 1 1/2" or 2" diameter access port at \$290
Certification	CE, UL/CSA (optional)

### Standard package

Part description	Quantity	Part image
AI 1.9 cu ft forced air oven	1 pc	
Stainless steel rack and brackets	2 sets	
110 volt 3-prong power cord	1 pc	
User's manual	1 pc	