

## GAS41 Polyclonal Antibody

Catalog number: 11050-1-AP

Size: 56 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

YEATS4; 4930573H17Rik,


B230215M10Rik, GAS41,

Glioma amplified sequence 41,

NUBI 1, NuB11, NuMA binding

protein 1, YAF9, YEATS domain

containing 4, YEATS4



human brain tissue were subjected to SDS PAGE followed by western blot with 11050-1-AP(YEATS4 antibody) at dilution of 1:300

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated GAS41 MW:	27 kDa
Observed GAS41 MW:	27-30 kDa
Positive WB detected in	Human brain tissue, A549 cells
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag1524
GenBank accession number:	BC000994
Gene ID (NCBI):	8089
Full name:	YEATS domain containing 4

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.