

HOXC4 Polyclonal Antibody

Catalog number: 14321-1-AP

Size: 48 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

HOXC4; HOX3; HOX3E; cp19



Jurkat cells were subjected to SDS PAGE followed by western blot with 14321-1-AP (HOXC4 Antibody) at dilution of 1:800

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Jurkat cells, HeLa cells
Calculated HOXC4 MW:	30kd
Observed HOXC4 MW:	30kd
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag5563
GenBank accession number:	BC050442
Gene ID (NCBI):	3221
Full name:	Homeobox C4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.