

## GABRG1 Polyclonal Antibody

Catalog number: 12871-1-AP

Size: 56 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

GABRG1; DKFZp686H2042,


GABA A receptor gamma 1,

GABAA receptor gamma 1,

GABRG1



Immunohistochemical of paraffin-embedded human gliomas using 12871-1-AP(GABRG1 antibody) at dilution of 1:100 (under 40x lens)



mouse eye tissue were subjected to SDS PAGE followed by western blot with 12871-1-AP(GABRG1 antibody) at dilution of 1:600

### Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated GABRG1 MW:	465aa, 54 kDa
Observed GABRG1 MW:	54 kDa
Positive WB detected in	Mouse eye tissue
Positive IHC detected in	Human gliomas tissue
Recommended dilution:	WB: 1:200-1:2000 IHC: 1:20-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3535
GenBank accession number:	BC031087
Gene ID (NCBI):	2565
Full name:	Gamma-aminobutyric acid (GABA) A receptor, gamma 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.