

## GCAP1 Polyclonal Antibody

Catalog number: 12944-1-AP

Size: 67 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

GUCA1A; C6orf131, COD3,

GCAP, GCAP 1, GCAP1,

Guanylate cyclase activator 1A,

GUCA, GUCA1, GUCA1A



Y79 cells were subjected to SDS PAGE followed by western blot with 12944-1-AP(GUCA1A antibody) at dilution of 1:200

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated GCAP1 MW:	23 kDa
Observed GCAP1 MW:	23-24kd
Positive WB detected in	Y79 cells, MCF7 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4017
GenBank accession number:	BC031663
Gene ID (NCBI):	2978
Full name:	Guanylate cyclase activator 1A (retina)

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.