

PGK2 Polyclonal Antibody

Catalog number: 13686-1-AP

Size: 48 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

PGK2; dJ417L20.2, PGK 2,

PGK2, PGKB, PGKPS,

phosphoglycerate kinase 2



mouse testis tissue were subjected to SDS PAGE followed by western blot with 13686-1-AP(PGK2 antibody) at dilution of 1:1200

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PGK2 MW:	417aa,45 kDa
Observed PGK2 MW:	45 kDa
Positive WB detected in	Mouse testis tissue, HeLa cells, HepG2 cells, human heart tissue, human testis tissue, mouse brain tissue
Recommended dilution:	WB: 1:1000-1:10000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4553
GenBank accession number:	BC038843
Gene ID (NCBI):	5232
Full name:	Phosphoglycerate kinase 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.