

JAM2 Polyclonal Antibody

Catalog number: 12972-1-AP

Size: 58 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

JAM2; C21orf43; CD322; JAM-B;

JAMB; PRO245; VE-JAM;

VEJAM



human placenta tissue were subjected to SDS PAGE followed by western blot with 12972-1-AP(JAM2 Antibody) at dilution of 1:1000

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human placenta tissue, human heart tissue
Calculated JAM2 MW:	298aa,33kd
Observed JAM2 MW:	33kd,43kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3431
GenBank accession number:	BC017779
Gene ID (NCBI):	58494
Full name:	Junctional adhesion molecule 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.