

MSRB3 Polyclonal Antibody

Catalog number: 14251-1-AP

Size: 55 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

MSRB3; DKFZp686C1178,

FLJ36866, MSRB3



mouse testis tissue were subjected to SDS PAGE followed by western blot with 14251-1-AP(MSRB3 Antibody) at dilution of 1:600

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated MSRB3 MW:	21 kDa
Observed MSRB3 MW:	26 kDa
Positive WB detected in	Mouse testis tissue
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC= Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag5524
GenBank accession number:	BC040053
Gene ID (NCBI):	253827
Full name:	Methionine sulfoxide reductase B3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.