

## MEIS3 Polyclonal Antibody

Catalog number: 12775-1-AP


Size: 61 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

MEIS3; DKFZp547H236; MRG2



mouse liver tissue were subjected to SDS PAGE followed by western blot with 12775-1-AP (MEIS3 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Mouse liver tissue
Calculated MEIS3 MW:	421aa,46kd
Observed MEIS3 MW:	46kd
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

### Immunogen information

Immunogen:	Ag3471
GenBank accession number:	BC025404
Gene ID (NCBI):	56917
Full name:	Meis homeobox 3

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.