

GABRA3 Polyclonal Antibody

Catalog number: 12708-1-AP


Size: 28 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

GABRA3; GABA A receptor
alpha 3, GABAA receptor alpha
3, GABRA3



human brain tissue were subjected to SDS PAGE followed by western blot with 12708-1-AP(GABRA3 antibody) at dilution of 1:100

Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated GABRA3 MW:	492aa,55 kDa
Observed GABRA3 MW:	55 kDa
Positive WB detected in	Human brain tissue, rat brain tissue
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3389
GenBank accession number:	BC028629
Gene ID (NCBI):	2556
Full name:	Gamma-aminobutyric acid (GABA) A receptor, alpha 3

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.