

alpha-SNAP Polyclonal Antibody

Catalog number: 10546-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

NAPA; alpha-SNAP, NAPA,

SNAP alpha, SNAPA

Background

SNAPs (soluble NSF-attachment proteins) may play an important role in the process of docking and fusion of vesicles to their target membranes, according to the 'SNARE hypothesis'. NAPA(Alpha- SNAP) was found in a wide range of tissues and act synergistically in intra-Golgi transport. This antibody recognizes a 33kd protein in SDS page.

Applications

Tested applications: ELISA

Species specificity: Human, Mouse, Rat; other species not tested.

Calculated alpha-SNAP MW: 33kd

Observed alpha-SNAP MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen: Ag0836

GenBank accession number: BC007432

Gene ID (NCBI): 8775

Full name: N-ethylmaleimide-sensitive factor attachment protein, alpha

Product information

Purification method: Antigen affinity purification

Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.