

PCDH1 Polyclonal Antibody

Catalog number: 13645-1-AP

Size: 29 µg/150 µl

Source: Rabbit

Isotype: IgG


Synonyms:

PCDH1; Cadherin like protein 1,

FLJ53887, PC42, PCDH1,

PCDH42, protocadherin 1,

Protocadherin 42



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 13645-1-AP (PCDH1 Antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PCDH1 MW:	1060aa, 115 kDa
Observed PCDH1 MW:	115 kDa
Positive WB detected in:	Fetal human brain tissue, human brain tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4591
GenBank accession number:	BC035812
Gene ID (NCBI):	5097
Full name:	Protocadherin 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.