



Model	TI-CS-T11	TI-CS-T1	JXG-213-F	JXG-513
Scan Type	Full Body	Full Body	Facial	Wrist
Setup Design	Tripod or Wall/Ceiling Mount	Tripod or handheld	Kiosk	Kiosk
Network Connection	Yes	No	Yes	Yes
Built-in Display	No	Yes	Yes	Yes
External Monitor Connection	Yes	Yes	No	No
Connect to PC	Yes	No	No	No
Alarm System	Yes	Yes	Yes	Yes
Power Type	Power Outlet	Battery or Power Outlet	Power Outlet	Power Outlet
Measurement Distance	up to 20FT	up to 15Ft	up to 3FT	up to 3CM
Max Person Scan	10 persons simultaneously	1 person at a time	1 person at a time	1 person at a time
Facial Recognition Capabilities	Yes	No	Yes	No
Built-in DB	Yes	No	Yes	Yes
Custom Temp. Range	Yes	Yes	Yes	Yes
Degree Error Range	-+ 0.3 C	-+ 0.5C	-+0.1 C	-+ 0.1C
API and Software	Yes	No	Yes	Yes
Video Formats	H.265, H.264, MJPEG	N/A	N/A	N/A
Storage	no SD card	Up to 8 GB SD Card	4GB	N/A
Networking	100M/1000M Ethernet, RJ45 interface Supports HTTP, TCP, RTSP, RTP, UDP, RTCP, Support ONVIF 28181 protocol	N/A	LAN×1, Wiegand Input×1, Wiegand Output×1, RS485×1	N/A
Product Includes	Software for data management - Need to provide own PC running Windows 10	Real time imaging Wireless (with rechargeable battery) or connect to outlet via power cable	LCD Display, Floorstand, cables, temperature sensor	LCD Display, Floorstand, cables, temperature sens
Price	\$15,900	\$12,900	\$1,599	\$1,399

**18-month warranty comes with all products*