

## PTPN9 Polyclonal Antibody

Catalog number: 12109-1-AP

Size: 40 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

PTPN9; MEG2, PTPase MEG2,

PTPMEG2, PTPN9



mouse liver tissue were subjected to SDS PAGE followed by western blot with 12109-1-AP( PTPN9 Antibody) at dilution of 1:1000

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Calculated PTPN9 MW:	593aa,68 kDa
Observed PTPN9 MW:	68 kDa
Positive WB detected in	Mouse liver tissue, A549 cells, mouse pancreas tissue, Transfected HEK-293 cells
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag2754
GenBank accession number:	BC010863
Gene ID (NCBI):	5780
Full name:	Protein tyrosine phosphatase, non-receptor type 9

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.