

HPCAL4 Polyclonal Antibody

Catalog number: 12470-1-AP


Size: 39 µg/150 µl

Source: Rabbit


Isotype: IgG

Synonyms:

HPCAL4; DKFZp761G122; HLP4



human brain tissue were subjected to SDS PAGE followed by western blot with 12470-1-AP(HPCAL4 Antibody) at dilution of 1:2000



Immunohistochemistry of paraffin-embedded human medulloblastoma using 12470-1-AP(HPCAL4 Antibody) at Dilution 1:100 (under 10x lens)

Applications

Tested applications:	ELISA, WB, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Positive WB detected in:	Human brain tissue
Calculated HPCAL4 MW:	191aa,22kd
Observed HPCAL4 MW:	22kd
Positive IHC detected in:	Human medulloblastoma
Recommended dilution:	WB: 1:1000-1:10000 IHC: 1:50-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag3153
GenBank accession number:	BC030827
Gene ID (NCBI):	51440
Full name:	Hippocalcin like 4

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.