

SRp20 PolyAb Antibody

Catalog number: 10916-1-AP

Size: 20 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

SFRS3; SRp20

Background

SFRS3, also named as SRSF3 and SRP20, belongs to the splicing factor SR family. It may be involved in RNA processing in relation with cellular proliferation and/or maturation. SFRS3 is one of the cognate factors. It is a proto-oncogene critical for cell proliferation and tumor induction and maintenance. The MW of SFRS3 is 20kd.

Applications

Tested applications: ELISA
Species specificity: Human, Mouse, Rat; other species not tested.
Calculated SRp20 MW: 20kd
Observed SRp20 MW:

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen: Ag1341
GenBank accession number: BC000914
Gene ID (NCBI): 6428
Full name: Splicing factor, arginine/serine-rich 3

Product information

Purification method: Antigen affinity purification
Storage: PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.