

# HyCare Blue Nitrile Examination Gloves



## Description:

Powder Free Nitrile Examination Gloves (M = 4.5+/- 0.2 gm). These blue powder-free nitrile examination gloves are ideal for a variety of applications. They provide superior fit and feel. This glove has exceptional tactile sensitivity and still offers excellent durability and puncture resistance.

## FEATURES

- Textured fingertips
- Strong chemical resistance
- High tactile comfort, yet strong
- Non-Sterile
- 1.5 AQL
- 510K Approved for medical use
- Approved for use with chemotherapy drugs
- Meet the requirements of ASTM D6319-10 and EN 455

## RECOMMENDED APPLICATIONS

- Medical Examinations
- Laboratory Research
- Food Handling
- Janitorial Work
- Light Duty Maintenance

## SPECIFICATIONS

MATERIAL:	100% Nitrile
POWDER-FREE:	Yes
COLOR:	Blue
GLOVE LENGTH:	240mm
FINGERTIP SURFACE:	Textured
FINGER THICKNESS:	4.35 mil
PALM THICKNESS:	3.15 mil
AVAILABLE SIZES:	S, M, L, XL
INSPECTION LEVEL (FFH):	AQL 1.5
NOT MADE FROM NATURAL RUBBER LATEX:	Yes

## PACKAGING

S-L: 100 gloves/box  
10 boxes/case

XL: 100 gloves/box  
10 boxes/case

## ORDERING INFORMATION

Thomas No.	Size	Qty.
20A00V269	S	100/Bx;1000/Cs
20A00V270	M	100/Bx;1000/Cs
20A00V271	L	100/Bx;1000/Cs
20A00V272	XL	100/Bx;1000/Cs

## DESIGN FEATURES:

Ambidextrous, straight fingers, finger textured and beaded cuff.

## QUALITY ASSURANCE:

Manufacturing process in compliance with ISO 9001 & ISO 13485 Quality Management System and ISO 14001 Environmental Management System.