

# Thermo Scientific Heraeus Megafuge 8 centrifuge



## Fits in. Stands out.

**NOW AVAILABLE**  
**APPLICATION-FOCUSED PACKAGES –**

*Order your centrifuge, rotor and accessories with 1 part number!*

**Fits in.** With its compact footprint, the Thermo Scientific™ Heraeus™ Megafuge™ 8 small benchtop centrifuge is sure to suit your space, and your processes and protocols, too.

**Stands-out** with the following innovative Thermo Scientific™ technologies:



- ▶ **Auto-Lock™ rotor exchange:**  
Secure push-button application versatility and cleaning convenience




- ▶ **ClickSeal™ biocontainment lids:**  
One-handed, certified<sup>1</sup> sample protection



- ▶ **Maximized capacity:**  
From 8 x 50 mL conical tubes to 24 x 5/7 mL blood tubes, and from microplates and microtubes, all in one versatile unit

<sup>1</sup> Biocontainment certification by Public Health England, Porton Down, UK.

# The Thermo Scientific Heraeus Megafuge 8 centrifuge fits in to your applications:




**RESEARCH APPLICATIONS**





- Cellular Biology
- Microbiology
- Genomics / Molecular Biology
- Proteomics
- Biochemistry
- Pharmaceutical Studies

**CLINICAL APPLICATIONS**




- Clinical Chemistry
- Clinical Microbiology
- Hematology
- Immunology
- Clinical Studies



## Stand-out Packages for Research Applications

Description	Thomas No.	Package Includes	
<b>Cell Culture:</b> up to 4 x 50 mL conical tubes	<b>1203P12</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>TX-150 swinging bucket rotor</li> <li>Round buckets (set of 4)</li> <li>Adapters for 50 mL conical or skirted tube (set of 4)</li> <li>Adapters for 15 mL conical tube (set of 4)</li> </ul>	
<b>High Capacity Cell Culture:</b> up to 8 x 50 mL conical tubes	<b>1203P13</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>TX-150 swinging bucket rotor</li> <li>50 mL conical buckets (unsealed, no adapter needed) (set of 4)</li> <li>Adapters for 15 mL conical tube (set of 4)</li> </ul>	
<b>High Capacity Cell Culture and Microtube:</b> up to 8x50 mL conical tubes and 30 x 2 mL microtubes	<b>1203P14</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>TX-150 swinging bucket rotor</li> <li>50 mL conical buckets (unsealed, no adapter needed) (set of 4)</li> <li>MicroClick 30 x 2 rotor</li> </ul>	
<b>Microplate:</b> up to 4 standard or 2 deepwell microplates	<b>1203P15</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>M10 swinging bucket rotor</li> <li>Unsealed buckets (set of 2)</li> </ul>	

## Stand-out Packages for Clinical Applications

Description	Cat. No.	Package Includes	
<b>High Capacity, Sealed:</b> up to 24 x 5/7 mL blood collection tubes	<b>1203P16</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>TX-150 swinging bucket rotor</li> <li>Round buckets (set of 4)</li> <li>ClickSeal biocontainment lids (set of 4)</li> <li>Adapters for 5/7 mL blood collection tube (set of 4)</li> <li>Adapters for 10 mL blood collection or 15 mL Corex™ / Kimble™ tube (set of 4)</li> </ul>	
<b>Medium Capacity:</b> up to 16 x 5/7 mL blood collection tubes	<b>1203P17</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>TX-100 clinical swinging bucket rotor with carriers</li> <li>Adapters for 5/7 mL BD Hemogard™ / BD Vacutainer™ tube (set of 4)</li> </ul>	
<b>Low Capacity, Sealed:</b> up to 8 x 5/7 mL blood collection tubes	<b>1203P18</b>	<ul style="list-style-type: none"> <li>Heraeus Megafuge 8 centrifuge, ventilated</li> <li>TX-100S clinical swinging bucket rotor with sealed carriers</li> <li>Adapters for 5/7 mL BD Hemogard™ / BD Vacutainer tube (set of 4)</li> </ul>	

Application-focused packages are available in the United States between January 1, 2014 and June 30, 2014 while supplies last. All offers subject to Thermo Fisher standard terms and conditions of sale. Not valid in connection with any other offer or discount. Void where prohibited by law or company policy. See your sales representative for additional information. Please note that the value of the specific offer(s) the customer may receive from the manufacturer under this program is a "discount or other reduction in price" to the customer under section 1128B(b)(3)(A) of the Social Security Act [42 U.S.C. 1320a-7b (b)(3)(a)]. Accordingly, the customer shall disclose this and any other discounts or other reduction in price received under this program under any state or federal program which provides cost or charge-based reimbursements to the customer for the products or services purchased under this program.

The Project was defined for JQC, Regulation 862.2050 – Clinical Centrifuges.

© 2014 Thermo Fisher Scientific Inc. All rights reserved. HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. Corex is a registered trademark of Corning Inc. Kimble is a registered trademark of Kimble Glass. BD Hemogard and BD Vacutainer are registered trademarks of BD Biosciences. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



**We Believe You Are Important**  
— How Can We Help?

800.345.2100 • 800.345.5232 fax  
PO Box 99 • Swedesboro, NJ 08085  
[www.thomassci.com](http://www.thomassci.com) • [value@thomassci.com](mailto:value@thomassci.com)



**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific