

Specification Sheet

PRODUCT: Hanks Balanced Salt Solution (HBSS)

CATALOG NO.: 472909 and
492909

LOT NO.: TBD

MANUFACTURE DATE TBD

EXPIRATION DATE: 1 Year from the date
of manufacture

STORAGE CONDITIONS: 4 - 40°C

Test Results

Test	Specification	Test Method	Result
Appearance	Clear/Colorless Liquid	Visual	TBD
pH	6.70 - 7.80	Potentiometric	TBD
Osmolality	270 - 305 mOsm	Freeze Point Depression	TBD
Protease	None Detected	Enzymatic	TBD
DNase	None Detected	Fluorometric	TBD
RNase	None Detected	Fluorometric	TBD
Peroxides Content	≤ 5ppb	Enzymatic	TBD
Cell Culture	Pass	Mouse Embryo Assay	TBD
Sterility	No Growth	USP <71>	TBD
Endotoxins	≤ 1.0 EU/mL	USP <85>	TBD

Not for direct human, animal or therapeutic use. For Research and Manufacturing Uses.

Manufacturing Site

Lifecycle Biotechnologies, 5858 Wright Drive, Loveland, Colorado 80538

Simple Lifecycle