

PHPT1 Polyclonal Antibody

Catalog number: 13112-1-AP

Size: 78 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:


PHPT1; bA216L13.10, CGI 202,

DKFZp564M173, HSPC141,

phosphohistidine phosphatase

1, PHP14, PHPT1, Protein janus

A homolog



MCF7 cells were subjected to SDS PAGE followed by western blot with 13112-1-AP (PHPT1 antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated PHPT1 MW:	125aa, 14 kDa
Observed PHPT1 MW:	14kd
Positive WB detected in	MCF7 cells, human kidney tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3778
GenBank accession number:	BC024648
Gene ID (NCBI):	29085
Full name:	Phosphohistidine phosphatase 1

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.