

HSF2 Polyclonal Antibody

Catalog number: 10706-1-AP

Size: 37 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

HSF2; MGC117376;

MGC156196; MGC75048



Immunohistochemistry of paraffin-embedded human prostate cancer using 10706-1-AP (HSF2 Antibody) at Dilution 1:100 (under 25x lens)

Applications

Tested applications:	ELISA, IHC
Species specificity:	Human, mouse, rat; other species not tested.
Calculated HSF2 MW:	79kd
Observed HSF2 MW:	
Positive IHC detected in:	Human prostate cancer
Recommended dilution:	IHC: 1:50-1:200

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation, FC = Flow cytometry

Immunogen information

Immunogen:	Ag1067
GenBank accession number:	BC005329
Gene ID (NCBI):	3298
Full name:	Heat shock transcription factor 2

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.