LRBFS20W1HC

Mediline





LIEBHERR

c Lister Us



Laboratory refrigerator

with spark-free interior

Total gross/net capacity I (cu.ft.) **Exterior dimensions** in mm (w/d/h)

in inches (w/d/h)

Interior dimensions in mm (w/d/h)

in inches (w/d/h)

Energy consumption in 365 days *

Ambient temperature in °C (°F)

Refrigerant

Sound rating **

Connection rating

Frequency/voltage

Cooling system

Defrost

Temperature range in °C (°F)

External cabinet finish

Door/lid material

Internal liner material

Type of control

Temperature display

Interior light

Adjustable shelves

Shelf material

Shelf load weight in kg (lbs)

Usable shelf area in mm (inches)

Handle

Lock

Door hinges

Insulation in mm (inches)

Gross/net weight in kg (lbs)

Accessories

Glass shelf

Additional lock barrels (up to 10)

LRBFS20W1HC

Mediline

567/533 (19.9/18.7) 747/743/1641 29 7/16"/29 5/16"/64 5/8" 600/560/1452 23 5/8"/22 1/16"/57 3/16" 380 kWh 10°C to 35°C (50°F to 95°F) R290 50 dB(A) 3.0 A

3.0 A 60 Hz/115V~ forced-air

automatic

1°C to 15°C (34°F to 59°F)

steel/white

SIEEI

commercial grade molded polystyrol

analog control external digital

n/a

5 glass

40 (88.2)

600/550 (23 5/8"/21 11/16") ergonomic slimline handle

recessed


right, reversible on site 32–52 (1 5/16"–2 1/16")

88/82 (194/180.8)

9293595

upon request

Unit Dimensions



Interior is free from sources of ignition and meets the requirements for Class I, Division 2 for storage of flammable substances in closed containers. Product complies with the applicable requirements of NFPA 45 and NFPA 99.

 $^{^\}star$ Measured unloaded at 26.6°C (80°F) ambient temperature with set point of 3.3°C (38°F).

 $^{^{\}star\star}$ Noise level – sound power re 1pW – according to IEC 60704-2-14.