

## SOX5 Polyclonal Antibody

Catalog number: 13216-1-AP

Size: 31 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

SOX5; L SOX5, SOX5,

Transcription factor SOX 5



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 13216-1-AP(SOX5 antibody) at dilution of 1:400

### Applications

Tested applications:	ELISA, WB
Cited applications:	WB
Species specificity:	Human,Mouse,Rat; other species not tested.
Cited species:	Mouse
Calculated SOX5 MW:	763aa,84 kDa
Observed SOX5 MW:	68-70 kDa
Positive WB detected in	SH-SY5Y cells, A549 cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag4092
GenBank accession number:	BC029220
Gene ID (NCBI):	6660
Full name:	SRY (sex determining region Y)-box 5

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.