

EDGE BIOLOGICALS, Inc.
598 North Second Street
Memphis, TN 38105
Phone: 901-523-0034
Fax: 901-527-3343

Certificate of Analysis

QUALITY ASSURANCE STATEMENT

Quality control testing has been performed on each lot of culture medium according to Document M22-A3 Approved Standard, "Quality Assurance for Commercially Prepared Microbiological Culture Media", published by the Clinical and Laboratory Standards Institute (CLSI).

PRODUCT: Sterile Saline 0.85% 10 mL		LOT #: 074	
PREPARED DATE: 03/15/2018		EXP. DATE: 03/15/2019	
pH: 7.27	APPEARANCE: Clear	METHOD of STERILIZATION: 121° C, 20 min	
Prepared from: DI Water	Lot #: N/A	Exp. Date: N/A	

This medium can be adjusted or supplemented according to the performance criteria required.

Quality Control

Sterility Test: Samples of this material were found to be sterile after three days incubation at 25± 2 degrees C and 35± 2 degrees C. No Growth observed at 72 hours in API Petri Film.

Release Date: 03/20/2018

Christopher Studaway
Quality Assurance Manager/Microbiologist

Christopher Studaway