

## Available in Standard or Hexagonal shapes

STERILE R

Manufactured and packaged in a clean room

### Sterile Weighing Boats

Ideal for applications where there is concern of cross contamination and 'clean' is a necessity in life science research, food, bio-pharmaceutical, environmental and molecular research environments

- Electron beam sterilized
- Lot number traceability for quality assurance
- Certificate of Sterility available upon request



Small



Small



Medium



Medium

Actual sizes

## Sterile Weighing Boats

Pack of 200/40 x 5pk

Shape	Regular Material Item No.	Antistatic Material Item No.	Color	Size	Volume mL	L x W x H in	L x W x H cm
Standard	—	120880	White <input type="checkbox"/>	Small	10	1.8 x 1.8 x .3	4.6 x 4.6 x .8
Standard	—	120710	White <input type="checkbox"/>	Medium	85	3.4 x 3.4 x .9	8.5 x 8.5 x 2.4
Hexagonal	—	120881	White <input type="checkbox"/>	Small	6	1.5 x 1 x .4	3.8 x 2.5 x .9
Hexagonal	—	120882	White <input type="checkbox"/>	Medium	46	2.5 x 2 x .7	6.4 x 5.1 x 1.9