

UGDH Polyclonal Antibody

Catalog number: 13151-1-AP

Size: 45 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:


UGDH; GDH, UDP Glc

dehydrogenase, UDP GlcDH,

UDP glucose 6 dehydrogenase,

UDP glucose dehydrogenase,

UDPGDH, UGD, UGDH



HepG2 cells were subjected to SDS PAGE followed by western blot with 13151-1-AP(UGDH antibody) at dilution of 1:500

Applications

Tested applications:	ELISA, WB
Cited applications:	IHC, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Cited species:	Human
Calculated UGDH MW:	494aa, 55 kDa
Observed UGDH MW:	55 kDa
Positive WB detected in	HepG2 cells, HeLa cells
Recommended dilution:	WB: 1:200-1:2000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag3767
GenBank accession number:	BC022781
Gene ID (NCBI):	7358
Full name:	UDP-glucose dehydrogenase

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.