

A35 NITRILE EXAMINATION GLOVES



FEATURES	SPECIFICATIONS	PACKAGING																								
<ul style="list-style-type: none"> • 4.3mil Thickness • Powder-Free • Beaded Cuff • Non-Sterile • Not Made From Natural Rubber Latex 	<table border="1"> <tr> <td>MATERIAL:</td> <td>Nitrile</td> </tr> <tr> <td>POWDER-FREE:</td> <td>Yes</td> </tr> <tr> <td>COLOR:</td> <td>Blue</td> </tr> <tr> <td>GLOVE LENGTH:</td> <td>245mm / 9.5 inches</td> </tr> <tr> <td>FINGERTIP SURFACE:</td> <td>Textured</td> </tr> <tr> <td>FINGER THICKNESS:</td> <td>0.11mm / 4.3mil</td> </tr> <tr> <td>PALM THICKNESS:</td> <td>0.08mm / 3.1mil</td> </tr> <tr> <td>AVAILABLE SIZES:</td> <td>XS, S, M, L, XL</td> </tr> <tr> <td>INSPECTION LEVEL (FFH):</td> <td>AQL 1.5</td> </tr> <tr> <td>NOT MADE FROM NATURAL RUBBER LATEX:</td> <td>Yes</td> </tr> </table>	MATERIAL:	Nitrile	POWDER-FREE:	Yes	COLOR:	Blue	GLOVE LENGTH:	245mm / 9.5 inches	FINGERTIP SURFACE:	Textured	FINGER THICKNESS:	0.11mm / 4.3mil	PALM THICKNESS:	0.08mm / 3.1mil	AVAILABLE SIZES:	XS, S, M, L, XL	INSPECTION LEVEL (FFH):	AQL 1.5	NOT MADE FROM NATURAL RUBBER LATEX:	Yes	<table border="1"> <tr> <td>XS-L: 100 gloves/box 10 boxes/case</td> </tr> <tr> <td>XL: 90 gloves/box 10 boxes/case</td> </tr> </table>	XS-L: 100 gloves/box 10 boxes/case	XL: 90 gloves/box 10 boxes/case		
MATERIAL:	Nitrile																									
POWDER-FREE:	Yes																									
COLOR:	Blue																									
GLOVE LENGTH:	245mm / 9.5 inches																									
FINGERTIP SURFACE:	Textured																									
FINGER THICKNESS:	0.11mm / 4.3mil																									
PALM THICKNESS:	0.08mm / 3.1mil																									
AVAILABLE SIZES:	XS, S, M, L, XL																									
INSPECTION LEVEL (FFH):	AQL 1.5																									
NOT MADE FROM NATURAL RUBBER LATEX:	Yes																									
XS-L: 100 gloves/box 10 boxes/case																										
XL: 90 gloves/box 10 boxes/case																										
RECOMMENDED APPLICATIONS		ORDERING INFORMATION																								
<ul style="list-style-type: none"> • Medical Examinations • Laboratory Research • Food Handling • Janitorial Work • Light Duty Maintenance 		<table border="1"> <thead> <tr> <th>Thomas No.</th> <th>Angenex No.</th> <th>Size</th> <th>Qty.</th> </tr> </thead> <tbody> <tr> <td>21A00A543</td> <td>A350</td> <td>XS</td> <td>100/Bx;1000/Cs</td> </tr> <tr> <td>21A00A544</td> <td>A351</td> <td>S</td> <td>100/Bx;1000/Cs</td> </tr> <tr> <td>21A00A545</td> <td>A352</td> <td>M</td> <td>100/Bx;1000/Cs</td> </tr> <tr> <td>21A00A546</td> <td>A353</td> <td>L</td> <td>100/Bx;1000/Cs</td> </tr> <tr> <td>21A00A547</td> <td>A354</td> <td>XL</td> <td>90/Bx;900/Cs</td> </tr> </tbody> </table>	Thomas No.	Angenex No.	Size	Qty.	21A00A543	A350	XS	100/Bx;1000/Cs	21A00A544	A351	S	100/Bx;1000/Cs	21A00A545	A352	M	100/Bx;1000/Cs	21A00A546	A353	L	100/Bx;1000/Cs	21A00A547	A354	XL	90/Bx;900/Cs
Thomas No.	Angenex No.	Size	Qty.																							
21A00A543	A350	XS	100/Bx;1000/Cs																							
21A00A544	A351	S	100/Bx;1000/Cs																							
21A00A545	A352	M	100/Bx;1000/Cs																							
21A00A546	A353	L	100/Bx;1000/Cs																							
21A00A547	A354	XL	90/Bx;900/Cs																							
REGULATORY COMPLIANCE																										
Cleared by FDA for medical use.																										
Component materials comply with FDA regulations for food contact (21 CFR 177-2600).																										