

CDH23 Polyclonal Antibody

Catalog number: 13496-1-AP

Size: 33 µg/150 µl

Source: Rabbit

Isotype: IgG

Synonyms:

CDH23; CAD23, Cadherin 23,

cadherin like 23, CDH23,

DFNB12, DKFZp434P2350,

FLJ00233, FLJ36499, KIAA1774,

KIAA1812, Otocadherin,

USH1D



A431 cells were subjected to SDS PAGE followed by western blot with 13496-1-AP (CDH23 antibody) at dilution of 1:300

Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated CDH23 MW:	3354aa, 369kd
Observed CDH23 MW:	59kd
Positive WB detected in	A431 cells, Y79 cells
Recommended dilution:	WB: 1:200-1:1000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

Immunogen information

Immunogen:	Ag4322
GenBank accession number:	BC032581
Gene ID (NCBI):	64072
Full name:	Cadherin-like 23

Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.02% sodium azide and 50% glycerol pH 7.3. Store at -20°C.