

Hei-CHILL 350 NR

Technical datasheet



Dimensions (W × H × D)	240 x 500 x 400 mm
Weight	31 kg
Display	LED
Working temperature range	-10°C...40°C
Display resolution	0.1°C
Interfaces	Alarm Output, RS 232
Resolution	N/A
Pump flow max.	9 L/min
Noise Level	57 dB(A)
Pump pressure max.	0.4 bar
Power Supply	115 V; 60 Hz
Filling Volume (Min)	4 L
Filling Volume (Max)	7 L
Power Consumption (Max)	0.5 kW
Power Plug	Cold Device Plug (IEC 60320, C14)
Protection Rating	IP 32
Temperature Stability	0.5 ± K
Refrigerant Stage 1	R-290 (GWP 3) A3; 0.041 kg; 0.0 t CO ₂ - eq
In/Outlet	8 mm
Safety Fittings	NFL

Resolution of Setting	0.1°C
Permissible Temperature Range During Storage	5°C...40°C
Permissible Temperature Range During Transport	-20°C...43°C
Protection Class	1