



KIMBLE® KIMAX® GRIFFIN BEAKERS, LOW FORM

BUY ONE GET ONE FREE MIX & MATCH FROM LIST

Designed for reliable everyday use, KIMBLE® KIMAX® low form beakers are ideal for routine mixing, heating, and general lab work. Made from durable borosilicate glass 3.3, they offer excellent heat and chemical resistance, easy pouring, clear graduations, and a practical marking area for everyday convenience.



KIMBLE® KIMAX® GRIFFIN BEAKERS

Thomas No.	Manufacturer No.	Capacity (mL)	OD (mm)	Height (mm)	Graduation Intervals (mL)	Graduation Range (mL)	Quantity
1530L26	14000-10	10	25	34	Not graduated	Not graduated	48
1530L46	14000-20	20	32	41	5	5 to 15	48
1530L56	14000-30	30	35	53	10	5 to 25	48
1530L66	14000-50	50	42	53	10	10 to 40	48
1530L76	14000-100	100	50	70	10	20 to 80	48
1530L86	14000-150	150	57	81	20	20 to 140	48
1530M11	14000-250	250	68	88	25	25 to 200	48
1530M21	14000-400	400	77	106	25	50 to 325	48
1530M31	14000-600	600	88	122	50	50 to 500	36
1530M41	14000-800	800	99	136	50	50 to 750	24
1530M51	14000-1000	1000	107	145	50	100 to 1000	24
1530M60	14000-1500	1500	120	165	100	200 to 1400	16
1530M71	14000-2000	2000	130	190	100	200 to 1800	8
1530M90	14000-4000	4000	161	250	250	500 to 3500	6

To claim your redemption (redeem by 12/31/2026)

1. Obtain your Proof of Purchase (dated between 05/01/26 - 12/31/26). Packing slip, purchase order and invoice are all considered valid for proof of purchase.
2. Visit www.dwkdeals.com
3. Attach your Proof of Purchase and fill-out the online form.
4. The value of freeware will be calculated off of DWK list price and units of measure found on www.dwk.com. Free items must be of equal or lesser value. Not for resale. Limited to two redemptions per promotion, per customer and/or per laboratory.

Offer valid for US and Canadian customers only. Availability for certain products may be limited by federal, state, provincial or local licensing requirements. Offer valid while supplies last. Not valid on prior purchases. DWK distributors and employees are prohibited from participating. Void in government accounts or where prohibited by law or company policy. Prices subject to change without notice. Trademarks contained herein are the property of their respective owners.

© 2026 DWK Life Sciences • The trademarks used are owned by DWK Life Sciences GmbH, DWK Life Sciences LLC or DWK Life Sciences Ltd.