

## VASH1 Polyclonal Antibody

Catalog number: 12730-1-AP

Size: 43 µg/150 µl


Source: Rabbit

Isotype: IgG

Synonyms:

VASH1; KIAA1036, VASH,

VASH1, vasohibin 1



human brain tissue were subjected to SDS PAGE followed by western blot with 12730-1-AP(VASH1 antibody) at dilution of 1:500

### Applications

Tested applications:	ELISA, WB
Species specificity:	Human, Mouse, Rat; other species not tested.
Calculated VASH1 MW:	365aa, 41 kDa
Observed VASH1 MW:	41 kDa
Positive WB detected in	Human brain tissue
Recommended dilution:	WB: 1:500-1:5000

Application key: WB = Western blotting, IHC = Immunohistochemistry, IF = Immunofluorescence, IP = Immunoprecipitation

### Immunogen information

Immunogen:	Ag3435
GenBank accession number:	BC009031
Gene ID (NCBI):	22846
Full name:	Vasohibin 1

### Product information

Purification method:	Antigen affinity purification
Storage:	PBS with 0.1% sodium azide and 50% glycerol pH 7.3. Store at -20°C.