


**HALYARD\* PUREZERO\*  
ULTRA VIOLET\***  
**Nitrile Exam Gloves**  
For Medical/Industrial/Lab Use

# TECHNICAL DATA SHEET

## Description

HALYARD\* **PUREZERO\*** ULTRA VIOLET\* (9.5") Nitrile Exam Gloves are designed for use in laboratories, research environments and clean areas with applications in pharmaceutical, medical device manufacturing, biotechnology and food processing/handling. HALYARD\* **PUREZERO\*** ULTRA VIOLET\* Nitrile Exam Gloves are tested for use against 12 chemotherapy drugs. The beaded cuff allows for ease of donning, and textured fingertips provide tactile agility.




## Physical Properties

AQL	1.0
Non-Sterile	✓
Ambidextrous	✓
Textured Fingertips	✓
Not Made With Natural Rubber Latex	✓
Powder-Free	✓
Tensile Strength <sup>1</sup>	20 MPa (Target)
Ultimate Elongation <sup>1</sup>	530%
Shelf Life	3 Years

## SMARTPULL\* Dispenser

Both **PUREZERO\*** ULTRA VIOLET\* gloves have a SMARTPULL\* Dispenser, which incorporates two separate openings on the box. The first, smaller opening is used when the box is full to control dispensing and lessen waste. When the box is half empty, the second, larger opening allows easier access to the gloves.



## Glove Dimensions

	X-SMALL	SMALL	MEDIUM	LARGE	X-LARGE
<b>ULTRA-VIOLET* Glove Length (inch/cm)</b>	9.5"/24	9.5"/24	9.5"/24	9.5"/24	9.5"/24
<b>Middle Finger Length (mm)</b>	69	77	78	85	88
<b>Width of Palm (mm)</b>	76	86	98	107	115
<b>Fingertip Thickness</b>	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)	0.09 mm (3.5 mil)
<b>Palm Thickness</b>	0.07 mm (2.8 mil)	0.07 mm (2.8 mil)	0.07 mm (2.8 mil)	0.07 mm (2.8 mil)	0.07 mm (2.8 mil)
<b>Cuff Thickness</b>	0.05 mm (2 mil)	0.05 mm (2 mil)	0.05 mm (2 mil)	0.05 mm (2 mil)	0.05 mm (2 mil)



**HALYARD\* PUREZERO\*  
ULTRA VIOLET\***  
Nitrile Exam Gloves  
For Medical/Industrial/Lab Use

# TECHNICAL DATA SHEET

## Ordering Information

Product Name	Thomas No.	Mfr. No.	Size	Box Qty.	Case Qty.
HALYARD* PUREZERO* ULTRA VIOLET* Nitrile Exam Gloves - 9.5"/24cm Length	22A00F347	LFSS11XS	XS	250	2500
	22A00F348	LFSS11SM	S	250	2500
	22A00F349	LFSS11MD	M	250	2500
	22A00F350	LFSS11LG	L	250	2500
	22A00F351	LFSS11XL	XL	230	2300

## Quality & Regulatory Standards

Compliant to these regulatory standards:

ISO 9001

ISO 13485: 2016

Compliance to REGULATION (EU) 2017/745

Compliant to these food handling regulatory standards:

FDA 21 CFR 177-2600

Commission Regulation (EU) No 10/2011

FDA 21 CFR part 820 accreditation

CE 2797 PPE Category III according to Regulation (EU) 2016/425 EEC

EN ISO 374-5:2016 Virus Protection

EN ISO 374-1:2016/Type C K-Low Chemical Protection

### US Standards

ASTM F739-12

ASTM D6978-05

### EU MDR Standards

EN 455-1:2020

EN 455-2:2015

EN 455-3:2015

EN 455-4:2009

### EU PPE

EN ISO 21420:2020

EN ISO 374-1:2016+A1:2018

EN ISO 374-2:2019

EN 16523-1:2015+A1:2018

EN ISO 374-4:2019

ISO-16604:2004



**CE 2797**  
(PPE Cat. III)

**ISO 374-5:2016**  
VIRUS

**UK CA**  
0086

**ISO 374-1:2016/Type C**  
K-Low Chemical



For additional information or samples,  
contact your local distributor or visit  
[www.purezerogloves.com](http://www.purezerogloves.com)

<sup>1</sup> Tested per ASTM D6319, EN 455-2

This fact sheet has been created using the most recent information. In the interest of continuous improvement, the characteristics of the products may change without prior notice.

\*Registered Trademark or Trademark of O&M Halyard or its affiliates. 2022. All rights reserved.  
008390 COPY-04125 June 2022



HALYARD\* is an Owens & Minor brand



ThomasSci.com  
833.544.SHIP (7447)  
CustomerService@thomassci.com

Connect With Us:

